



SummitHEALTHCARE

SUMMIT CARE EXCHANGE

CLINICAL EXCELLENCE THROUGH CONTINUITY OF CARE

Meeting transition of care requirements not only leads to improved patient outcomes, but is also a federally mandated Meaningful Use stage two objective. The continuity of care document implements a standard for the data exchange to ensure the information is accepted and usable by all participating systems.

The Summit Care Exchange solution facilitates Continuity of Care Document (CCD) exchange among hospitals, physician practices, and HIEs to promote interoperability. It allows for the seamless exchange of data to support the transitions of patient care. Summit Care Exchange supports the three CCD integration models; Discharge, Consumer/Supplier and Direct Messaging using the Summit Healthcare HISP.

PRODUCT DETAILS

Benefits

- Support your transition of care meaningful use objective
- Support better clinical decision making for improved patient outcomes
- Enable clinicians to send critical electronic medical information to other providers without loss of meaning
- Configure your CCD interfaces for easy connection to multiple other entities
- One centralized place for monitoring and alerting of CCD interface traffic
- Support Direct Messaging including HISP offering

Features

- Configure **inbound/outbound CCD Integration connections** with multiple entities
- **Send and consume CCDs** with another entity for support of the MU transition of care requirement
- **Share all file and report types** with your providers
- Fully integrated with the **Summit Dashboard** which allows centralized viewing for monitoring and alerting
- Integrated Summit Healthcare HISP for Direct Messaging
- Query HIE and RLS Repositories for Patients CCDs
- Ingest CCD's with the added flexibility to filter or remove inactive, duplicate or old data

www.summit-healthcare.com

35 Braintree Hill Park, Suite 303, Braintree, MA 02184

Phone: (781) 519-4840

Email: info@summit-healthcare.com

SUMMIT CARE EXCHANGE

SUMMIT CARE EXCHANGE DASHBOARD



* Summit Care Exchange with Summit Healthcare Direct supports XDS.b, ITI, SOAP+XDR/XDM, SMIME/SMTP, Direct+XDR/XDM, MTOM, PIX.

SUMMIT HEALTHCARE PRODUCT SUITE



Summit Express Connect

A comprehensive interface engine that provides dashboard monitoring & alerting. Seamlessly integrate your EHR with other vendor systems.



Summit Provider Exchange

An enhanced bi-directional physician office integration solution, which simplifies patient registration, inbound orders and outbound results between your hospital EHR and ambulatory systems.



Summit Care Exchange

Streamline and support your transition of care meaningful use objective by supporting your CCD data exchange. Provides better support, and more informed clinical decision making.



Summit Scripting Toolkit

A workflow and task automation tool that can also be used when HL7 methods to deliver data are not an option. Has the ability to handle everything from the simplest script to the most complex interface.



Summit InSync

The market leading solution for dictionary/table management through extraction, analysis and synchronization. Keep your data accurate between test and live, across multiple platforms and domains.



Summit Downtime Reporting System

Be prepared for any hospital EHR, Ancillary system or Network downtime with a cost effective, simple, yet robust technology. Ensure your staff has access to critical patient data 24/7.

PRODUCT REQUIREMENTS

Operating System	Microsoft Windows XP, Microsoft Windows Vista, or Microsoft Windows 7
Minimum Memory	2 GB RAM
Minimum Hard Disk Space	50 MB
Minimum System Processor	2 gigahertz (GHz) or faster 32-bit (x86) or 64-bit (x64) processor
Required Software	.NET Framework 3.5 or higher, Microsoft SQL Server 2005 or above, Standard or Enterprise Edition, Summit Express Connect
Minimum Resolution	1024x768